

#

---

```
# This is a bidsmap YAML file with the key-value mappings for the
different BIDS
# datatypes (anat, func, dwi, etc). The datatype attributes are the
keys that map
# onto the BIDS labels. The bidsmap data-structure should be 5
levels deep:
# dict > dict > list > dict > dict
#
# NB:
# 1) Edit the bidsmap file to your needs before feeding it to
bidscoiner.py
# 2) (Institute) users may create their own bidsmap_[template].yaml
or
# bidsmap_[sample].yaml file
#
# For more information, see:
# https://github.com/Donders-Institute/bidscoin
# https://docs.ansible.com/ansible/latest/reference_appendices/
YAMLSyntax.html
#
```

---

#### Options:

```
  bidscoin:
    version: 3.0.8 # BIDScoin version (should
correspond with the version in ../bidscoin/version.txt)
    bidsignore: extra_data/ # Semicolon-separated list of
entries that are added to the .bidsignore file (for more info, see
BIDS specifications), e.g. extra_data;/pet;/myfile.txt;yourfile.csv
    dcm2niix: # See dcm2niix -h and https://
www.nitrc.org/plugins/mwiki/index.php/dcm2nii:MainPage#General_Usage
for more info
    path: /Users/ah2816/opt/anaconda3/envs/Python/bin/dcm2niix; #
Command to set the path to dcm2niix (note the semi-colon), e.g.
module add dcm2niix/1.0.20180622; or PATH=/opt/dcm2niix/bin:$PATH;
or /opt/dcm2niix/bin/ or "C:\Program Files\dcm2niix\" (note the
quotes to deal with the whitespace)
    args: -b y -z y -i n # Argument string that is passed
to dcm2niix. Tip: SPM users may want to use '-z n' (which produces
unzipped nifti's, see dcm2niix -h for more information)
```

#### PlugIns: []

##### DICOM:

```
  subject: <<SourceFilePath>> # <<SourceFilePath>> extracts the
subject label from the source directory during bidscoiner runtime. A
DICOM attribute can also be used as subject-label instead, e.g.
<<PatientID>>
  session: <<SourceFilePath>> # <<SourceFilePath>> extracts the
```

session label from the source directory during bidscoiner runtime. A DICOM attribute can also be used as session-label instead, e.g. <StudyID>

```
anat:      # ----- All anatomical runs
-----
func:      # ----- All functional runs
-----
perf:      # ----- All perfusion runs
-----
```

```
- provenance:
  attributes: &perf_dicomattr
    Modality:
    ProtocolName:
    SeriesDescription:
    ImageType:
    SequenceName:
    SequenceVariant:
    ScanningSequence:
    MRAcquisitionType:
    SliceThickness:
    FlipAngle:
    EchoNumbers:
    EchoTime:
    RepetitionTime:
    PhaseEncodingDirection:
  bids:
    acq: <ProtocolName>
    rec: <SeriesDescription>
    dir:
    run: <<1>>
    suffix: asl
- provenance:
  attributes: *perf_dicomattr
  bids:
    <<:
```

# See: schema/datatypes/

```
perf.yaml
  acq: <ProtocolName>
  rec: <SeriesDescription>
  dir:
  run: <<1>>
  suffix: asl
  suffix: m0scan
  acq: <ProtocolName>
  rec: <SeriesDescription>
  dir:
  run: <<1>>
- provenance:
  attributes: *perf_dicomattr
  bids:
    acq: <ProtocolName>
    rec: <SeriesDescription>
    dir:
```

```

    run: <<1>>
    suffix: aslcontext
- provenance:
  attributes: *perf_dicomattr
  bids:
    acq: <ProtocolName>
    rec: <SeriesDescription>
    run: <<1>>
    suffix: asllabeling

dwi:      # ----- All diffusion runs
-----
fmap:     # ----- All fieldmap runs
-----
eeg:      # ----- All EEG data
-----
- provenance:
  attributes:
    Modality: EEG
    ProtocolName:
    SeriesDescription:
  bids:                                     # See: schema/datatypes/
eeg.yaml
  acq: <ProtocolName>
  task: <SeriesDescription>
  run: <<1>>
  suffix: eeg

ieeg:     # ----- All iEEG data
-----
- provenance:
  attributes:
    Modality: EEG
    ProtocolName:
    SeriesDescription:
  bids:                                     # See: schema/datatypes/
ieeg.yaml
  acq: <ProtocolName>
  task: <SeriesDescription>
  run: <<1>>
  suffix: ieeg

beh:      # ----- All behavioural data
-----
pet:      # ----- All PET runs
-----
extra_data: # ----- All extra data
-----
leave_out: # ----- Data that will be left out
-----
PAR:
  subject: <<SourceFilePath>> # <<SourceFilePath>> extracts the
subject label from the source directory during bidscoiner runtime. A
DICOM attribute can also be used as subject-label instead, e.g.

```

<PatientID>  
session: <<SourceFilePath>> # <<SourceFilePath>> extracts the session label from the source directory during bidscoiner runtime. A DICOM attribute can also be used as session-label instead, e.g.  
<StudyID>

anat: # ----- All anatomical runs

-----  
- provenance: /Volumes/toshiba/ICCAM/DICOMS/3/sub-3002/ses-1/  
sub-3002\_T1w.PAR

attributes:

series\_type: Image MRSERIES  
exam\_name: MRH050  
protocol\_name: T1W\_IR 1150 SENSE  
tech: T1TFE  
diffusion: 0  
fov: '[240. 144. 191.25]'  
scan\_resolution: '[256 205]'  
scan\_mode: 3D  
max\_slices: 160  
flow\_compensation: 0  
max\_echoes: 1  
repetition\_time: '[8.4]'  
prep\_direction: Right-Left

bids:

acq: ''  
ce:  
rec:  
run: ''  
part: ['', mag, phase, real, imag, 0]  
suffix: T1w

- provenance: /Volumes/toshiba/ICCAM/DICOMS/3/sub-3002/ses-1/  
sub-3002\_T2w.PAR

attributes:

series\_type: Image MRSERIES  
exam\_name: MRH050  
protocol\_name: Anna\_PD SENSE  
tech: TSE  
diffusion: 0  
fov: '[256. 153.65 192. ]'  
scan\_resolution: '[256 195]'  
scan\_mode: MS  
max\_slices: 40  
flow\_compensation: 0  
max\_echoes: 1  
repetition\_time: '[6980.]'  
prep\_direction: Right-Left

bids:

suffix: PDT2  
acq: ''  
ce:  
rec:  
run: ''  
part: ['', mag, phase, real, imag, 0]

```
func:          # ----- All functional runs
-----
- provenance: /Volumes/toshiba/ICCAM/DICOMS/3/sub-3002/ses-1/
sub-3002_evo_run-1.PAR
attributes:
  series_type: Image  MRSERIES
  exam_name: MRH050
  protocol_name: EVOC run1 2E SENSE
  tech: FEEPI
  diffusion: 0
  fov: '[225. 99. 225.]'
  scan_resolution: '[64 61]'
  scan_mode: MS
  max_slices: 33
  flow_compensation: 0
  max_echoes: 2
  repetition_time: '[2000.001]'
  prep_direction: Anterior-Posterior
bids:
  task: evo
  acq: ''
  ce:
  dir:
  rec:
  run: <<1>>
  echo:
  part: ['', mag, phase, real, imag, 0]
  suffix: bold
- provenance: /Volumes/toshiba/ICCAM/DICOMS/3/sub-3002/ses-1/
sub-3002_evo_run-2.PAR
attributes:
  series_type: Image  MRSERIES
  exam_name: MRH050
  protocol_name: EVOC run2 2E SENSE
  tech: FEEPI
  diffusion: 0
  fov: '[225. 102. 225.]'
  scan_resolution: '[64 61]'
  scan_mode: MS
  max_slices: 34
  flow_compensation: 0
  max_echoes: 2
  repetition_time: '[1999.999]'
  prep_direction: Anterior-Posterior
bids:
  task: evo
  acq: ''
  ce:
  dir:
  rec:
  run: <<2>>
  echo:
  part: ['', mag, phase, real, imag, 0]
  suffix: bold
```

- provenance: /Volumes/toshiba/ICCAM/DICOMS/3/sub-3002/ses-1/  
sub-3002\_gng\_run-1.PAR

attributes:

series\_type: Image MRSERIES  
exam\_name: MRH050  
protocol\_name: GO\_NO\_GO run1 2E SENSE  
tech: FEEPI  
diffusion: 0  
fov: '[225. 99. 225.]'  
scan\_resolution: '[64 61]'  
scan\_mode: MS  
max\_slices: 33  
flow\_compensation: 0  
max\_echoes: 2  
repetition\_time: '[2000.001]'  
prep\_direction: Anterior-Posterior

bids:

task: gng  
acq: ''  
ce:  
dir:  
rec:  
run: <<1>>  
echo:  
part: ['', mag, phase, real, imag, 0]  
suffix: bold

- provenance: /Volumes/toshiba/ICCAM/DICOMS/3/sub-3002/ses-1/  
sub-3002\_gng\_run-2.PAR

attributes:

series\_type: Image MRSERIES  
exam\_name: MRH050  
protocol\_name: GO\_NO\_GO run2 2E SENSE  
tech: FEEPI  
diffusion: 0  
fov: '[225. 99. 225.]'  
scan\_resolution: '[64 61]'  
scan\_mode: MS  
max\_slices: 33  
flow\_compensation: 0  
max\_echoes: 2  
repetition\_time: '[2000.001]'  
prep\_direction: Anterior-Posterior

bids:

task: gng  
acq: ''  
ce:  
dir:  
rec:  
run: <<2>>  
echo:  
part: ['', mag, phase, real, imag, 0]  
suffix: bold

- provenance: /Volumes/toshiba/ICCAM/DICOMS/3/sub-3002/ses-1/  
sub-3002\_mid\_run-1.PAR

```
attributes:
  series_type: Image    MRSERIES
  exam_name: MRH050
  protocol_name: MID run1 2E SENSE
  tech: FEEPI
  diffusion: 0
  fov: '[225. 99. 225.]'
  scan_resolution: '[64 61]'
  scan_mode: MS
  max_slices: 33
  flow_compensation: 0
  max_echoes: 2
  repetition_time: '[2000.001]'
  prep_direction: Anterior-Posterior
bids:
  task: mid
  acq: ''
  ce:
  dir:
  rec:
  run: <<1>>
  echo:
  part: ['', mag, phase, real, imag, 0]
  suffix: bold
- provenance: /Volumes/toshiba/ICCAM/DICOMS/3/sub-3002/ses-1/
sub-3002_mid_run-2.PAR
attributes:
  series_type: Image    MRSERIES
  exam_name: MRH050
  protocol_name: MID run2 2E SENSE
  tech: FEEPI
  diffusion: 0
  fov: '[225. 99. 225.]'
  scan_resolution: '[64 61]'
  scan_mode: MS
  max_slices: 33
  flow_compensation: 0
  max_echoes: 2
  repetition_time: '[2000.001]'
  prep_direction: Anterior-Posterior
bids:
  task: mid
  acq: ''
  ce:
  dir:
  rec:
  run: <<2>>
  echo:
  part: ['', mag, phase, real, imag, 0]
  suffix: bold
- provenance: /Volumes/toshiba/ICCAM/DICOMS/3/sub-3002/ses-1/
sub-3002_rest_run-1.PAR
attributes:
  series_type: Image    MRSERIES
```

```

exam_name: MRH050
protocol_name: RESTING_STATE SENSE
tech: FEEPI
diffusion: 0
fov: '[225. 102. 225.]'
scan_resolution: '[64 61]'
scan_mode: MS
max_slices: 34
flow_compensation: 0
max_echoes: 2
repetition_time: '[1999.999]'
prep_direction: Anterior-Posterior
 bids:
  task: rest
  acq: ''
  ce:
  dir:
  rec:
  run: <<1>>
  echo:
  part: ['', mag, phase, real, imag, 0]
  suffix: bold
perf:      # ----- All perfusion runs
-----
- provenance:
  attributes: &perf_parattr
  series_type:
  exam_name:
  protocol_name:
  tech:
  diffusion:
  fov:
  scan_resolution:
  scan_mode:
  max_slices:
  flow_compensation:
  max_echoes:
  repetition_time:
  prep_direction:
  bids:
    acq: <protocol_name>
    rec: <tech>
    dir:
    run: <<1>>
    suffix: asl
- provenance:
  attributes: *perf_parattr
  bids:
    <<:
# See: schema/datatypes/
perf.yaml
  acq: <protocol_name>
  rec: <tech>
  dir:

```



```

        run: <<1>>
        suffix: asl
    suffix: m0scan
    acq: <protocol_name>
    rec: <tech>
    dir:
    run: <<1>>
- provenance:
  attributes: *perf_parattr
  bids:
    acq: <protocol_name>
    rec: <tech>
    dir:
    run: <<1>>
    suffix: aslcontext
- provenance:
  attributes: *perf_parattr
  bids:
    acq: <protocol_name>
    rec: <tech>
    run: <<1>>
    suffix: asllabeling

dwi:          # ----- All diffusion runs
-----
fmap:         # ----- All fieldmap runs
-----
- provenance: /Volumes/toshiba/ICCAM/DICOMS/3/sub-3002/ses-1/
sub-3002_fmap.PAR
  attributes:
    series_type: Image    MRSERIES
    exam_name: MRH050
    protocol_name: WIP Anna_B0_mapping CLEAR
    tech: T1FFE
    diffusion: 0
    fov: '[225.      98.988 225.   ]'
    scan_resolution: '[64 64]'
    scan_mode: MS
    max_slices: 33
    flow_compensation: 0
    max_echoes: 1
    repetition_time: '[599.]'
    prep_direction: Anterior-Posterior
  bids:
    acq: ''
    run: ''
    suffix: fieldmap
    IntendedFor: ''
eeg:          # ----- All EEG data
-----

ieeg:        # ----- All iEEG data
-----

```

```

meg:          # ----- All MEG data
-----

beh:          # ----- All behavioural data
-----

pet:          # ----- All PET runs
-----

extra_data:   # ----- All extra data
-----

leave_out:    # ----- Data that will be left out
-----
- provenance: /Volumes/toshiba/ICCAM/DICOMS/3/sub-3002/ses-1/
sub-3002_DTI.PAR
  attributes:
    series_type: Image  MRSERIES
    exam_name: MRH050
    protocol_name: WIP DTI_Laurence_64dirs SENSE
    tech: DwiSE
    diffusion: 1
    fov: '[240. 133. 240.]'
    scan_resolution: '[128 124]'
    scan_mode: MS
    max_slices: 70
    flow_compensation: 0
    max_echoes: 1
    repetition_time: '[9299.997]'
    prep_direction: Anterior-Posterior
  bids:
    acq: <protocol_name><exam_name>
    rec: DwiSE
    suffix: ImageMRSERIES
- provenance: /Volumes/toshiba/ICCAM/DICOMS/3/sub-3002/ses-1/
sub-3002_localizer.PAR
  attributes:
    series_type: Image  MRSERIES
    exam_name: MRH050
    protocol_name: Survey
    tech: T1TFE
    diffusion: 0
    fov: '[250. 250. 50.]'
    scan_resolution: '[256 128]'
    scan_mode: MS
    max_slices: 9
    flow_compensation: 0
    max_echoes: 1
    repetition_time: '[11.125]'
    prep_direction: Anterior-Posterior
  bids:
    acq: Survey
    inv:
    rec: T1TFE

```

```
ce:
task:
echo:
dir:
run: <<1>>
mod:
part:
suffix:
```

```
P7:
#
```

```
-----
# P*.7 key-value heuristics (GE fields that are mapped to the BIDS
labels)
#
```

```
Nifti:
#
```

```
-----
# Nifti key-value heuristics (Nifti fields that are mapped to the
BIDS labels)
#
```

```
FileSystem:
#
```

```
-----
# File system key-value heuristics (these file- and folder names will
be mapped
# to the BIDS labels; Special substitutions can be performed using
python's
# Format Specification Mini-Language)
#
```